



An Occlusion Can Be Irreversibly Damaged by Using Occlusal Contact Indicators.

Philip L Millstein DMD, MS*

*Department of Restorative Dentistry, Harvard School of Dental Medicine, Boston, Massachusetts, USA.

Corresponding author: Philip L Millstein DMD, MS, Department of Restorative Dentistry, Harvard School of Dental Medicine, Boston, Massachusetts, USA, E-mail: pmillsteindmd@gmail.com

Abstract

The improper usage of occlusal contact indicators can result in destructive tooth reductions. New tooth additions are not possible. The positive and negative effects of indicator markings and their interpretations must be addressed.

Keywords: Indicators, Contact, Force, Intensity.

Dentistry was formally established in 1867. Vulcanite rubber was used as the base material for dentures. Crude porcelain teeth were inserted for function and appearance. Function was minimal. Manual typewriters were introduced. The denture occlusion was monitored using typewriter ribbon and carbon paper (1). Commercial marketing of paper and ribbon followed. Not much has changed. Electronic devices may be employed to record occlusal contact forces however marking indicators are still required. Indicated accuracy is unknown (2). Twenty eight gauge indicator wax is used to measure contact (3). A black

wax pencil is used to mark exposed areas of teeth that perforate the green wax. Wax markings are not precise and they are often incomplete. There are no standards. Dentists use millions of indicators. They are disposable. The cost is minimal to the dentist compared to the profits reaped by the manufacturer. No clinician or dental testing organization questions indicator accuracy or reliability (4). All indicators make contact but the results differ [Figure 1]. Marking systems bleed ink upon contact. The electronic systems measure force [Figure 2].



Figure 1. All contact areas differ. Different markings are made by indicators. Size, shape, and intensity vary.

The pressed indicator wax attempts to display contact intensity. Knowledgeable clinicians know what they use and use what they know. They know their indicators. Some clinicians choose force indicating systems that do not mark? [Figure 3]. Is a force measuring system necessary or can contacts be recorded in simple fashion with a contact indicator? How much unnecessary tooth structure is removed in an attempt to balance a bite on a wafer using force measuring wafers? An alternative means is to make visual comparisons of indicated


tooth contact using a novel contact indicator (5). This requires taking a silicone impression of an occlusal contact. The contact surface does not distort. An indicator can damage a contact marking upon impact.

Contact areas can be highlighted using the Easy Scan with image analysis. The total procedure takes a minute. The exact size and location of the contact areas is shown in red. Contact areas are not conventional [Figure 4]. Intensity of contact is readable.

Received: October 20, 2026; **Revised:** October 23, 2026; **Accepted:** October 25, 2026

Citation: Millstein PL. (2026) An Occlusion Can Be Irreversibly Damaged by Using Occlusal Contact Indicators. J Oral Health Dent Res, 5(3): 1-3.

Copyright: ©2026 Millstein PL. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.



Available Digital Occlusal Technology Systems...

- T-Scan 10 Novus** Tekscan Inc, USA – *40 years of System Evolution.*
 Inception date 1984 – *256 Relative Occlusal Force Levels captured in 0.003 secs/frame*
- Occlusense I** Bausch GMBH, Germany –
5 years since inception - 256 Relative Occlusal Force Levels captured in 0.056 secs/frame
Marks teeth similar to Articulating Paper
- Accura I** Dmetec, South Korea –
5 years since inception – 256 Relative Occlusal Force Levels in 0.01 secs/frame
- IOS Scanners (Many)** – Optically scan occlusal morphology with all teeth out of contact
No Relative Occlusal Force Levels; no Timing except ModJaw

Misreading contact markings and force renderings in daily dental practice is common. Removing essential tooth structure occurs. The damage goes unnoticed. The clinician remains unaware. Occlusal equilibration is a common dental procedure. Insurance pays. Manufacturers sell a disposable item (6). Electronic systems with wafers are expensive and may not provide an improvement over existing methods. Without a testing standard there is no way to compare and evaluate measuring systems. Every occlusal contact is unique (7). No two are the same. Marking indicators tend to make contacts appear conventional. Contacts appear uniform (8). Schools teach ideal occlusal contact anatomy. Contact recognition is a strong factor in

choosing an indicator. Clinical success is dependent on marking interpretations. Who is to judge a true contact when there are no standards? A viable occlusion depends upon strong occlusal contact (9). Removing contacts threatens occlusal stability. How many implant reconstructions are the product of occlusal instability and a failed occlusion. Class room instruction on the choice and application of occlusal contact media could be instituted in the dental curriculum. Contact indicating is a significant part of clinical dentistry. It should be included under the discipline of preventive dentistry.




Figure 4. The impressed contact area is printed in shades of red. The contacts are not abbreviated, they are true. The areas appears in unconventional shapes. No two contact areas are alike. Contacts are contact areas.

REFERENCES

1. Millstein PL. An Evaluation of Occlusal Contact Marking Indicators: A Descriptive Quantitative Method. Quintessence International 1983;14:813
2. Occlusense Bausch GMBH, Germany
3. Occlusal Indicator Wax. Kerr Manufacturing Company. Brea, CA kerrdental.com
4. American Dental Association. Materials Testing and Specifications. Somerville MA
5. Millstein PL., Sabrosa CE., Florencio S., Geber K. Consistency and Repeatability of Digitized Occlusal Records IJRRD 3(2):24-29,2020; No.58152
6. Millstein PL. Occlusal Contact Marking Indicators. They All Differ. J Oral Health Dent Res, 5(3): 2026
7. Millstein PL., Sabrosa CE., and Wai Yung. Occlusal Damage Can Be Recorded. Do We Repair It? J Oral Health Dent Res, 5(3): 2026
8. L and G Steier The Come-Up Technique. Tufts School of Dental Medicine. Boston, MA . USA
9. Millstein PL., Wai Yung and Sabrosa CE. Posterior Bite Collapse. J Oral Health Dent Res, 5(2): 2026